

# SAFE® aspen fine

PRODUCT DATA SHEET

Release date: February 2025

## Definition

### ASPEN WOOD GRANULATE

Made from selected premium raw material (trunk wood)

## Product Purpose

Bedding material

## Directions for Use

### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

### SHELF-LIFE

36 months from date of production

### AUTOCLAVABLE



SAFE® aspen fine

Picture indicative only

## Particle Size Distribution

Average particle size/particle range		2 mm - 3 mm	
Sieve		> 2.5 mm	> 2 mm
		> 1.25 mm	> 0.71 mm
Retained	15 - 45 %	50 - 85 %	90 - 100 %
			97 - 100%

## Chemical and Physical Properties

SAFE® aspen fine	
Color	yellow
Structure	cubic granulate
Ash (850 °C, 4 h)	~ 1 %
Bulk density	~ 195 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

## Product Presentation

Standard Packaging: 10 kg autoclavable polypropylene bags.

Big bags and bulk trucks are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.

## Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany



Plant Germany



Pictures indicative only