

## SAFE® BK 8-15

PRODUCT DATA SHEET  
Release date: February 2025

### Definition

#### SPRUCE WOOD GRANULATE

Made from selected premium raw material (trunk wood)

### Product Purpose

Bedding material

### Directions for Use

#### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

#### SHELF-LIFE

36 months from date of production

#### AUTOCLAVABLE

### Particle Size Distribution

Average particle size/particle range	2.5 mm - 3.5 mm			
Sieve	> 2.5 mm	> 2 mm	> 1.25 mm	> 0.71 mm
Retained	25 - 50 %	60 - 85 %	min. 90 %	min. 97 %



SAFE® BK 8-15

Picture indicative only

### Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany



Plant Germany

### Chemical and Physical Properties

#### SAFE® BK 8-15

Color	yellow
Structure	cubic granulate
Ash (850 °C, 4 h)	~ 0.5 %
Bulk density	~ 1.95 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

### Product Presentation

Standard Packaging: 10 kg autoclavable polypropylene bags.

Big bags and bulk trucks are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.



Pictures indicative only