

SAFE® MK 2000

PRODUCT DATA SHEET

Release date: February 2025

Definition

CORNCOB GRANULATE

Made from selected premium raw material. The raw material used is non-genetically modified corncobs in accordance with EU regulations. (EC) No.1829/2003 and No.1830/2003.

Product Purpose

Bedding material

Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

SHELF-LIFE

36 months from date of production

AUTOCCLAVABLE

Particle Size Distribution

Average particle size/particle range	1.7 mm - 2.2 mm			
Sieve	> 2.5 mm	> 2 mm	> 1.6 mm	> 0.71 mm
Retained	max. 0.5 %	max. 30 %	min. 50 %	min. 99 %

Chemical and Physical Properties

SAFE® MK 2000

Color	light brown
Structure	granulate
Ash (850 °C, 4 h)	~ 1 %
Bulk density	~ 550 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

Product Presentation

Standard Packaging: 20 kg autoclavable polypropylene bags.

Big bags and bulk trucks are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.



SAFE® MK 2000

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany