

SAFE® performance^{plus}

PRODUCT DATA SHEET

Release date: February 2025

Definition

CELLULOSE PELLETS AND CELLULOSE FLAKES

Made from selected premium raw material

Product Purpose

Bedding and Enrichment material

Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

SHELF-LIFE

36 months from date of production

AUTOCLAVABLE

Particle Size Distribution

Average particle size/particle range 6 mm - 20 mm

Chemical and Physical Properties

SAFE® performance^{plus}

Color	light grey and white
Structure	pellets and flakes
Bulk density	~ 150 g/l
pH-value (10 % suspension)	5.5 - 7.5
100 % natural product	

Product Presentation

Standard Packaging: 35 l autoclavable paper bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.

SAFE® performance^{plus}

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany