

## SAFE<sup>®</sup> performance<sup>plus</sup>

PRODUCT DATA SHEET  
Release date: February 2025

### Definition

#### CELLULOSE PELLETS AND CELLULOSE FLAKES

Made from selected premium raw material

### Product Purpose

Bedding and Enrichment material

### Directions for Use

#### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

#### SHELF-LIFE

36 months from date of production

#### AUTOCLAVABLE

### Particle Size Distribution

Average particle size/particle range 6 mm - 20 mm

### Chemical and Physical Properties

SAFE<sup>®</sup> performance<sup>plus</sup>

Color	light grey and white
Structure	pellets and flakes
Bulk density	~ 150 g/l
pH-value (10 % suspension)	5.5 - 7.5
100 % natural product	

### Product Presentation

Standard Packaging: 35 l autoclavable paper bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.



SAFE<sup>®</sup> performance<sup>plus</sup>

Picture indicative only

### Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany