

# SAFE® pure

PRODUCT DATA SHEET

Release date: February 2025

## Definition

### CELLULOSE FLAKES

Made from selected premium raw material

## Product Purpose

Bedding material

## Directions for Use

### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

### SHELF-LIFE

36 months from date of production

### AUTOCCLAVABLE

## Chemical and Physical Properties

SAFE® pure

|                            |           |
|----------------------------|-----------|
| Color                      | white     |
| Structure                  | flakes    |
| Ash (850 °C, 4 h)          | max 2.5 % |
| Bulk density               | ~ 130 g/l |
| pH-value (10 % suspension) | 5.5 - 7.5 |
| 100 % natural product      |           |



SAFE® pure

Picture indicative only

## Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany

## Product Presentation

Standard Packaging: 7 kg autoclavable polypropylene bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.



Picture indicative only