

SAFE® 3-4 S

PRODUCT DATA SHEET

Release date: July 2019

Definition

SPRUCE WOOD FIBRILLATED FIBERS

Made from selected premium raw material (trunk wood)

Product Purpose

Bedding material

Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

SHELF-LIFE

36 months from date of production

Particle Size Distribution

Average particle size/particle range	2.0 mm - 4.0 mm		
Sieve	> 3.55 mm	> 2 mm	> 1.25 mm
Retained	max. 5 %	max. 57 %	min. 80 %



SAFE® 3-4 S

Picture indicative only

Certificates



Chemical and Physical Properties

SAFE® 3-4 S

Color	beige
Structure	fibrillated fibers
Ash (850° C, 4 h)	~ 0.5 %
Bulk density	~ 85 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

Product Presentation

Standard Packaging: 12.5 kg autoclavable polypropylene bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Analysis certificate of contaminants provided upon request.



Picture indicative only

Produced in Germany