

SAFE® gnawing small

PRODUCT DATA SHEET
Release date: September 2026

Definition

SELECTED WOOD

Product Purpose

Enrichment material

USE

Nesting	Gnawing	Activity
	x	x



Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from light

SHELF-LIFE

36 months from date of production

AUTOCLAVABLE

Product Properties

SAFE® gnawing small

Color	light brown
Structure	block
Dimensions	100 mm x 20 mm x 5 mm

SAFE® gnawing small

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany



Plant Germany

Product Presentation

Standard Packaging: ca. 500 units per carton

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.