

SAFE® comfort natural

PRODUCT DATA SHEET
Release date: February 2025

Definition

CELLULOSE FLAKES

Made from selected premium raw material

Product Purpose

Bedding and Enrichment material

Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

SHELF-LIFE

36 months from date of production

Particle Size Distribution

Average particle size/particle range	10 mm - 25 mm		
Sieve	> 10 mm	> 5 mm	> 2 mm
Retained	max. 35 %	min. 75 %	min. 95 %



SAFE® comfort natural

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany

Chemical and Physical Properties

SAFE® comfort natural

Color	light brown to grey
Structure	flakes
Bulk density	~ 55 g/l
pH-value (10 % suspension)	5.5 - 7.5
100 % natural product	

Product Presentation

Standard Packaging: 60 l non autoclavable PE-bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.