

# SAFE<sup>®</sup> compact crinklets natural

PRODUCT DATA SHEET  
Release date: February 2025

## Definition

SELECTED PREMIUM PAPER, COMPACTED

## Product Purpose

Enrichment material

### USE

Nesting	Gnawing	Activity
x		x



## Directions for Use

### STORAGE

Store in a clean and dry place with constant temperature, protected from light

### SHELF-LIFE

36 months from date of production

### AUTOCLAVABLE

## Product Properties

SAFE<sup>®</sup> compact crinklets natural

Color	brown
Structure	compacted crinklets
Dimensions	Ø ca. 50 mm; ca. 10 mm (H)

## Product Presentation

Standard Packaging: 700 units per carton

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.

SAFE<sup>®</sup> compact crinklets natural

Picture indicative only

## Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany