

SAFE[®] compact nesting small

PRODUCT DATA SHEET
Release date: February 2025

Definition

SELECTED SOFTWOOD SHAVINGS, COMPACTED

Product Purpose

Enrichment material

USE

Nesting	Gnawing	Activity
x		x



Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from light

SHELF-LIFE

36 months from date of production

AUTOCLAVABLE

Product Properties

SAFE[®] compact nesting small

Color	light brown
Structure	compacted shavings
Dimensions	Ø ca. 50 mm; ca. 10 mm (H)

Product Presentation

Standard Packaging: 500 units per carton

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.

SAFE[®] compact nesting small

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany



Plant Germany