

SAFE® block beech XL

PRODUCT DATA SHEET

Release date: February 2025

Definition


BEECH WOOD

Product Purpose

Enrichment material

USE

Nesting	Gnawing	Activity
	x	x



Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from light

SHELF-LIFE

36 months from date of production

AUTOCLAVABLE

Product Properties

SAFE® block beech XL

Color	light brown
Structure	block
Dimensions	96 mm x 40 mm x 40 mm

SAFE® block beech XL

Picture indicative only

Certificates



Plant Germany
Plant France



Plant Germany



Plant Germany



Plant Germany

Product Presentation

Standard Packaging: ca. 96 units per carton.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request