

# SAFE® comfort white

PRODUCT DATA SHEET

Release date: February 2025

## Definition

### CELLULOSE FLAKES

Made from selected premium raw material

## Product Purpose

Bedding and Enrichment material

## Directions for Use

### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

### SHELF-LIFE

36 months from date of production

## Particle Size Distribution

Average particle size/particle range		10 mm - 25 mm	
Sieve	> 10 mm	> 5 mm	> 2 mm
Retained	max. 35 %	min. 75 %	min. 95 %



SAFE® comfort white

Picture indicative only

## Certificates



Plant Germany  
Plant France



Plant Germany



Plant Germany

## Chemical and Physical Properties

SAFE® comfort white

Color	white
Structure	flakes
Bulk density	~ 55 g/l
pH-value (10 % suspension)	5.5 - 7.5
100 % natural product	

## Product Presentation

Standard Packaging: 50 l non autoclavable PE-bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.