

## SAFE® hut

PRODUCT DATA SHEET  
Release date: September 2024

### Definition

CARDBOARD HOUSE

### Product Purpose

Enrichment material

#### USE

Nesting	Gnawing	Activity
x		x



### Directions for Use

#### STORAGE

Store in a clean and dry place with constant temperature, protected from light

#### SHELF-LIFE

36 months from date of production

#### AUTOCLAVABLE

### Product Properties

SAFE® hut

Color	white
Structure	house
Dimensions	140 x 90 x 60 mm

SAFE® hut is manufactured from virgin pulp white paper.

### Product Presentation

Standard Packaging: 480 units per carton.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Quarterly analysis certificate of contaminants provided upon request.

SAFE® hut

Picture indicative only

### Certificates

